



Crustacean hyperglycemic hormones Antibody

Product Code	CSB-PA318620ZA01CDS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14944
Immunogen	Recombinant Carcinus maenas Crustacean hyperglycemic hormones protein
Raised In	Rabbit
Species Reactivity	Carcinus maenas
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Carcinus maenas (Common shore crab) (Green crab)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Crustacean hyperglycemic hormones,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.