



dxr Antibody

Product Code	CSB-PA332636ZA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45568
Immunogen	Recombinant Escherichia coli dxr protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Escherichia coli (strain K12)
Lead Time	Made-to-order (12-14 weeks)
Target Names	dxr
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.