



Odorant binding protein OBP28 Antibody

Product Code	CSB-PA31438ZA01DZM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B0XC79
Immunogen	Recombinant Culex quinquefasciatus Odorant binding protein OBP28 protein
Raised In	Rabbit
Species Reactivity	Culex quinquefasciatus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Culex quinquefasciatus (Southern house mosquito) (Culex pungens)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Odorant binding protein OBP28
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.