



DEFB106A Antibody

Product Code	CSB-PA711382ZA01EQV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5IAB3
Immunogen	Recombinant Pan troglodytes DEFB106A protein
Raised In	Rabbit
Species Reactivity	Pan troglodytes
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Pan troglodytes (Chimpanzee)
Lead Time	Made-to-order (12-14 weeks)
Target Names	DEFB106A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.