



hmfB Antibody

Product Code	CSB-PA322515ZA01MFF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19267
Immunogen	Recombinant Methanothermus fervidus hmfB protein
Raised In	Rabbit
Species Reactivity	Methanothermus fervidus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Methanothermus fervidus
Lead Time	Made-to-order (12-14 weeks)
Target Names	hmfB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.