



hmw3 Antibody

Product Code	CSB-PA690351ZA01MLW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q50360
Immunogen	Recombinant Mycoplasma pneumoniae hmw3 protein
Raised In	Rabbit
Species Reactivity	Mycoplasma pneumoniae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Mycoplasma pneumoniae (strain ATCC 29342 / M129)
Lead Time	Made-to-order (12-14 weeks)
Target Names	hmw3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.