



rev Antibody

Product Code	CSB-PA748326ZA01FYK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q75006
Immunogen	Recombinant Human immunodeficiency virus type 1 group M subtype C rev protein
Raised In	Rabbit
Species Reactivity	Human immunodeficiency virus type 1 group M subtype C
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Human immunodeficiency virus type 1 group M subtype C (isolate ETH2220) (HIV-1)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	rev