



ompC Antibody

Product Code	CSB-PA669669ZA01KBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q48473
Immunogen	Recombinant Klebsiella pneumoniae ompC protein
Raised In	Rabbit
Species Reactivity	Klebsiella pneumoniae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Klebsiella pneumoniae
Lead Time	Made-to-order (12-14 weeks)
Gene Names	ompC