




Major mite allergen Der p 23 Antibody

Product Code	CSB-PA31708ZA01DJH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	L7N6F8
Immunogen	Recombinant Dermatophagoides pteronyssinus Major mite allergen Der p 23 (22-90aa)
Species Reactivity	Dermatophagoides pteronyssinus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Dermatophagoides pteronyssinus
Target Names	Major mite allergen Der p 23

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Major mite allergen Der p 23 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 20 kDa

Observed band size: 23 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.