



pstP Antibody

Product Code	CSB-PA301778ZA01GKB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WHW4
Immunogen	Recombinant Mycobacterium tuberculosis pstP protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Lead Time	Made-to-order (12-14 weeks)
Target Names	pstP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.