



Genome polyprotein, partial Antibody

Product Code	CSB-PA31148ZA01GLA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F6KTB0
Immunogen	Recombinant Enterovirus A71 Genome polyprotein,partial protein
Raised In	Rabbit
Species Reactivity	Enterovirus A71
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Enterovirus A71
Lead Time	Made-to-order (12-14 weeks)
Target Names	Genome polyprotein,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.