



mxs Antibody

Product Code	CSB-PA812928ZA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8JH68
Immunogen	Recombinant Danio rerio mxs protein
Raised In	Rabbit
Species Reactivity	Danio rerio
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Lead Time	Made-to-order (12-14 weeks)
Target Names	mxs
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.