



hmfA Antibody

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Product Code | CSB-PA343755ZA01MFF |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P48781 |
| Immunogen | Recombinant Methanothermus fervidus hmfA protein |
| Raised In | Rabbit |
| Species Reactivity | Methanothermus fervidus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugate |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Immunogen Species | Methanothermus fervidus |
| Lead Time | Made-to-order (12-14 weeks) |
| Target Names | hmfA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |