



chp Antibody

Product Code	CSB-PA746696ZA01SKV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GFB3
Immunogen	Recombinant Staphylococcus aureus chp protein
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Staphylococcus aureus (strain MRSA252)
Lead Time	Made-to-order (12-14 weeks)
Target Names	chp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.