



Envelope protein p54 Antibody

Product Code	CSB-PA717644ZA01AEC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q65272
Immunogen	Recombinant African swine fever virus Envelope protein p54 protein
Raised In	Rabbit
Species Reactivity	African swine fever virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	African swine fever virus (isolate Pig/Portugal/Lis 60/1960) (ASFV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Envelope protein p54,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.