



ompA Antibody

Product Code	CSB-PA332912ZA01DSA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23421
Immunogen	Recombinant Chlamydia trachomatis ompA protein
Raised In	Rabbit
Species Reactivity	Chlamydia trachomatis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Chlamydia trachomatis
Lead Time	Made-to-order (12-14 weeks)
Target Names	ompA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.