



epd Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA367173ZA01AUH |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | A0KGD2 |
| Immunogen | Recombinant Aeromonas hydrophila subsp. hydrophila epd protein |
| Raised In | Rabbit |
| Species Reactivity | Aeromonas hydrophila subsp. hydrophila |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugate |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Immunogen Species | Aeromonas hydrophila subsp. hydrophila (strain ATCC 7966 / DSM 30187 / JCM 1027 / KCTC 2358 / NCIMB 9240) |
| Lead Time | Made-to-order (12-14 weeks) |
| Target Names | epd |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |