



INS Antibody

Product Code	CSB-PA011742ZA01HO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01310
Immunogen	Recombinant Equus caballus INS protein
Raised In	Rabbit
Species Reactivity	Equus caballus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Equus caballus (Horse)
Lead Time	Made-to-order (12-14 weeks)
Target Names	INS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.