



aroE Antibody

Product Code	CSB-PA683186ZA01SAAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5HNV1
Immunogen	Recombinant Staphylococcus epidermidis aroE protein
Raised In	Rabbit
Species Reactivity	Staphylococcus epidermidis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Staphylococcus epidermidis (strain ATCC 35984 / RP62A)
Lead Time	Made-to-order (12-14 weeks)
Target Names	aroE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.