



fbpA Antibody

Product Code	CSB-PA374958ZA01MOJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1KQD8
Immunogen	Recombinant Mycobacterium bovis fbpA protein
Raised In	Rabbit
Species Reactivity	Mycobacterium bovis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Mycobacterium bovis (strain BCG / Pasteur 1173P2)
Lead Time	Made-to-order (12-14 weeks)
Target Names	fbpA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.