



Major capsid protein VP1 Antibody

Product Code	CSB-PA402662ZA01WBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5HBD5
Immunogen	Recombinant WU polyomavirus Major capsid protein VP1 protein
Raised In	Rabbit
Species Reactivity	WU polyomavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	WU polyomavirus (WUPyV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Major capsid protein VP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.