



mutL Antibody

Product Code	CSB-PA309097ZA01TNN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P96082
Immunogen	Recombinant Thermus aquaticus mutL protein
Raised In	Rabbit
Species Reactivity	Thermus aquaticus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Thermus aquaticus
Lead Time	Made-to-order (12-14 weeks)
Gene Names	mutL