



crtX Antibody

Product Code	CSB-PA325132ZA01ERA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21686
Immunogen	Recombinant Pantoea ananas crtX protein
Raised In	Rabbit
Species Reactivity	Pantoea ananas
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Pantoea ananas (Erwinia uredovora)
Lead Time	Made-to-order (12-14 weeks)
Target Names	crtX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.