



groEL Antibody

Product Code	CSB-PA625985ZA01HAAH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1CVE5
Immunogen	Recombinant Helicobacter pylori groEL protein
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Helicobacter pylori (strain HPAG1)
Lead Time	Made-to-order (12-14 weeks)
Target Names	groEL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.