



acrR Antibody

Product Code	CSB-PA32608ZA01SGI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3Z4T5
Immunogen	Recombinant Shigella sonnei acrR protein
Raised In	Rabbit
Species Reactivity	Shigella sonnei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Shigella sonnei (strain Ss046)
Lead Time	Made-to-order (12-14 weeks)
Target Names	acrR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.