



Chi3l1 Antibody

Product Code	CSB-PA893408ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9WTV1
Immunogen	Recombinant Rattus norvegicus Chi3l1 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Chi3l1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.