



Bm1_01445 Antibody

Product Code	CSB-PA427056ZA01BWV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8NFF0
Immunogen	Recombinant Brugia malayi Bm1_01445 protein
Raised In	Rabbit
Species Reactivity	Brugia malayi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Brugia malayi (Filarial nematode worm)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Bm1_01445
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.