



acpS Antibody

Product Code	CSB-PA352442ZA01FVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63471
Immunogen	Recombinant Streptococcus agalactiae serotype III acpS protein
Raised In	Rabbit
Species Reactivity	Streptococcus agalactiae serotype III
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Streptococcus agalactiae serotype III (strain NEM316)
Lead Time	Made-to-order (12-14 weeks)
Target Names	acpS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.