



inaX Antibody

Product Code	CSB-PA325961ZA01XAJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18127
Immunogen	Recombinant Xanthomonas campestris pv. translucens inaX protein
Raised In	Rabbit
Species Reactivity	Xanthomonas campestris pv. translucens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Xanthomonas campestris pv. translucens
Lead Time	Made-to-order (12-14 weeks)
Target Names	inaX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.