



yopH Antibody

Product Code	CSB-PA325344ZA01YAAQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15273
Immunogen	Recombinant Yersinia enterocolitica yopH protein
Raised In	Rabbit
Species Reactivity	Yersinia enterocolitica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Yersinia enterocolitica
Lead Time	Made-to-order (12-14 weeks)
Target Names	yopH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.