



S4 Antibody

Product Code	CSB-PA365971ZA01RCH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03527
Immunogen	Recombinant Reovirus type 3 S4 protein
Raised In	Rabbit
Species Reactivity	Reovirus type 3
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Reovirus type 3 (strain Dearing) (T3D) (Mammalian orthoreovirus 3)
Lead Time	Made-to-order (12-14 weeks)
Target Names	S4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.