



sapA Antibody

Product Code	CSB-PA330493ZA01CYZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35827
Immunogen	Recombinant Campylobacter fetus sapA protein
Raised In	Rabbit
Species Reactivity	Campylobacter fetus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Campylobacter fetus
Lead Time	Made-to-order (12-14 weeks)
Gene Names	sapA