



sox1 Antibody

Product Code	CSB-PA022417ZA01XBF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A2TED3
Immunogen	Recombinant Xenopus tropicalis sox1 protein
Raised In	Rabbit
Species Reactivity	Xenopus tropicalis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Xenopus tropicalis (Western clawed frog) (Silurana tropicalis)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	sox1