



L3 Antibody

Product Code	CSB-PA334192ZA01HIG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36849
Immunogen	Recombinant Human adenovirus B serotype 3 L3 protein
Raised In	Rabbit
Species Reactivity	Human adenovirus B serotype 3
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Human adenovirus B serotype 3 (HAdV-3) (Human adenovirus 3)
Lead Time	Made-to-order (12-14 weeks)
Target Names	L3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.