



gB Antibody

Product Code	CSB-PA361840ZA01HWV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06473
Immunogen	Recombinant Human cytomegalovirus gB protein
Raised In	Rabbit
Species Reactivity	Human cytomegalovirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (Human herpesvirus 5)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	gB