



ophA1 Antibody

Product Code	CSB-PA327204ZA01BXR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33164
Immunogen	Recombinant Burkholderia cepacia ophA1 protein
Raised In	Rabbit
Species Reactivity	Burkholderia cepacia
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Burkholderia cepacia (Pseudomonas cepacia)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ophA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.