



elavl4 Antibody

Product Code	CSB-PA007595ZA01XBF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4QNI8
Immunogen	Recombinant Xenopus tropicalis elavl4 protein
Raised In	Rabbit
Species Reactivity	Xenopus tropicalis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Xenopus tropicalis (Western clawed frog) (Silurana tropicalis)
Lead Time	Made-to-order (12-14 weeks)
Target Names	elavl4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.