



GRA3 Antibody

Product Code	CSB-PA467284ZA01TOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B6KEU8
Immunogen	Recombinant Toxoplasma gondii GRA3 protein
Raised In	Rabbit
Species Reactivity	Toxoplasma gondii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Toxoplasma gondii
Lead Time	Made-to-order (12-14 weeks)
Target Names	GRA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.