



# N Antibody

Product Code	CSB-PA349523ZA01HQE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59595
Immunogen	Recombinant Human SARS coronavirus N protein
Raised In	Rabbit
Species Reactivity	Human SARS coronavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Human SARS coronavirus (SARS-CoV) (Severe acute respiratory syndrome coronavirus)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	N