



resilin Antibody

Product Code	CSB-PA892693ZA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9V7U0
Immunogen	Recombinant Drosophila melanogaster resilin protein
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Drosophila melanogaster (Fruit fly)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	resilin