



U100 Antibody

Product Code	CSB-PA865550ZA01HKA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9QJ11
Immunogen	Recombinant Human herpesvirus 6B U100 protein
Raised In	Rabbit
Species Reactivity	Human herpesvirus 6B
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Human herpesvirus 6B (strain Z29) (HHV-6 variant B) (Human B lymphotropic virus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	U100
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.