



U69 Antibody

Product Code	CSB-PA329474ZA01HJZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24446
Immunogen	Recombinant Human herpesvirus 6A U69 protein
Raised In	Rabbit
Species Reactivity	Human herpesvirus 6A
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Human herpesvirus 6A (strain Uganda-1102) (HHV-6 variant A) (Human B lymphotropic virus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	U69
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.