



U96/U97/U98/U99/U100 Antibody

Product Code	CSB-PA731651ZA01HJY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q69489
Immunogen	Recombinant Human herpesvirus 6A U96/U97/U98/U99/U100 protein
Raised In	Rabbit
Species Reactivity	Human herpesvirus 6A
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Human herpesvirus 6A (strain GS) (HHV-6 variant A) (Human B lymphotropic virus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	U96/U97/U98/U99/U100
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.