



Minor capsid protein VP2 Antibody

Product Code	CSB-PA333002ZA01MKK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24596
Immunogen	Recombinant Murine polyomavirus Minor capsid protein VP2 protein
Raised In	Rabbit
Species Reactivity	Murine polyomavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Murine polyomavirus (strain Kilham) (MPyV) (Murine pneumotropic virus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Minor capsid protein VP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.