



mip Antibody

Product Code	CSB-PA346473ZA01LDJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53605
Immunogen	Recombinant Legionella longbeachae mip protein
Raised In	Rabbit
Species Reactivity	Legionella longbeachae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Legionella longbeachae
Lead Time	Made-to-order (12-14 weeks)
Target Names	mip
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.