



tdnL Antibody

Product Code	CSB-PA842386ZA01FFZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93JW0
Immunogen	Recombinant Pseudomonas putida tdnL protein
Raised In	Rabbit
Species Reactivity	Pseudomonas putida
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Pseudomonas putida (Arthrobacter siderocapsulatus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	tdnL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.