



gyrB Antibody

Product Code	CSB-PA358232ZA01SWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A2I4
Immunogen	Recombinant Salmonella typhi gyrB protein
Raised In	Rabbit
Species Reactivity	Salmonella typhi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Salmonella typhi
Lead Time	Made-to-order (12-14 weeks)
Gene Names	gyrB