



4hbD Antibody

Product Code	CSB-PA327911ZA01CWX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38945
Immunogen	Recombinant Clostridium kluuyveri 4hbD protein
Raised In	Rabbit
Species Reactivity	Clostridium kluuyveri
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Clostridium kluuyveri (strain ATCC 8527 / DSM 555 / NCIMB 10680)
Lead Time	Made-to-order (12-14 weeks)
Target Names	4hbD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.