




L1R Antibody

Product Code	CSB-PA324949ZA01VAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20540
Immunogen	Recombinant Vaccinia virus Protein L1(L1R),partial (2-183aa)
Species Reactivity	Vaccinia virus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vaccinia virus
Target Names	L1R

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: L1R antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 26.8 kDa

Observed band size: 29 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.