



Heat shock 70 kDa Antibody

Product Code	CSB-PA362326ZA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08106
Immunogen	Recombinant Gallus gallus Heat shock 70 kDa protein
Raised In	Rabbit
Species Reactivity	Gallus gallus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Gallus gallus (Chicken)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Heat shock 70 kDa protein,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.